

ANTI A, ANTI B, ANTI AB, ANTI D BLOOD GROUPING REAGENT



DETECT THE PRESENCE OR ABSENCE OF A, B OR RHESUS ANTIGENS ON THE SURFACE OF HUMAN RED BLOOD CELLS

The Blood Grouping reagents are used for the identification of blood types. The test procedure is based on the agglutination principle, where red cells possessing the typing antigen agglutinate in the presence of the corresponding antibody in the testing reagent indicating the presence of the tested antigen. No agglutination indicates the absence of the tested antigen.

Atlas Medical ABO reagents are prepared from In-Vitro culture supernatants of hybridized immunoglobulin-secreting mouse cell lines. The reagents are formulated and optimized for use in tube and slide methods. Atlas Medical provides high quality blood grouping reagents that are accurate, easy to use, competitively priced, and conveniently packed in different sizes and options.

3

MINUTES
RESULT



EASY
TO USE



FAST RESULT

MDA APPROVED



Certificate No.: IVDD6168923-139572

PRODUCT INFORMATION

ITEM	QUANTITY	STORAGE & STABILITY
ATLAS Anti A Monoclonal Reagent	10ml/ vial	2°C – 8°C stability: until expiry
ATLAS Anti B Monoclonal Reagent		
ATLAS Anti AB Monoclonal Reagent		
ATLAS Anti D IgG/IgM Blend Reagent		

BLOOD GROUPING REAGENTS

- Anti A monoclonal reagent
- Anti B monoclonal reagent
- Anti AB monoclonal reagent
- Anti D IgG/IgM Blend reagent



RESULT INTERPRETATION

+ POSITIVE
If no agglutination is observed.

- NEGATIVE
If agglutination appears.

Anti A monoclonal reagent	Anti B monoclonal reagent	Anti AB monoclonal reagent	Anti D IgG/ IgM blend reagent	ABO Group
+	-	+	+	A+
+	-	+	-	A-
-	+	+	+	B+
-	+	+	-	B-
+	+	+	+	AB+
+	+	+	-	AB-
-	-	-	+	O+
-	-	-	-	O-

ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY
8.02.00.0.0010	ATLAS Anti A Monoclonal Reagent (Titer: 1 /512)	1 vial/ kit
8.02.01.0.0010	ATLAS Anti B Monoclonal Reagent (Titer: 1 /512)	1 vial/ kit
8.02.02.0.0010	ATLAS Anti AB Monoclonal Reagent (Titer: 1 /512)	1 vial/ kit
8.02.03.0.0010	ATLAS Anti D IgG/IgM Blend Reagent (Titer: 1 /128)	1 vial/ kit

DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

Address: B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.

Tel: +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my

