# MAS<sup>®</sup> Alcohol/Ammonia

## LIQUID ASSAYED ALCOHOL/AMMONIA CONTROL

IVD

### INTENDED USE

**MAS<sup>®</sup> Alcohol/Ammonia** is intended for use in the clinical laboratory as a quantitative control for monitoring test procedures used to assay human serum and/or plasma specimens for alcohol and/or ammonia. Include **Alcohol/Ammonia** with patient serum or plasma specimens when assaying for alcohol and/or ammonia. Assay values are provided for specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

### PRODUCT DESCRIPTION

**Alcohol/Ammonia** is a liquid stable control material prepared from human serum. Analyte levels are adjusted with various pure chemicals. Preservatives and stabilizers are added to maintain product integrity.

**DANGER: Alcohol/Ammonia** control contains ≤0.1% human serum albumin (HSA), and ≤0.001% 2-Methyl-2H-isothiazol-3-one  
H317 - May cause allergic skin reaction.  
H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

**Alcohol/Ammonia** is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009. The packaging of this product contains dry natural rubber.

### STORAGE AND STABILITY

Store **Alcohol/Ammonia** at 2-8°C. Unopened vials are stable until the expiration date on the label. Once opened, vials of control are stable for 30 days when stored tightly capped at 2-8°C. **Do not freeze.**

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

### CONTROL RANGES

The assigned ranges for these controls are based upon replicate assays of representative samples of the product by participating laboratories in accordance with established protocol. All values have been assigned with instruments and instrument manufacturer's reagents available at the time of assay. Subsequent instrument or reagent modifications may invalidate these assigned ranges.

Expected values may vary slightly with different reagent and/or methodologies used. Refer to the included table for values obtained for specific systems. Values listed are specific for this lot of control only. Good laboratory practice suggests that each laboratory establish its own parameters.

### INSTRUCTIONS FOR USE

Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C. Assay controls in accordance with the reagent manufacturer's instructions for unknown specimens.

Dropper tips are available to help minimize exposure of control samples to air. DO NOT use a syringe needle to withdraw sample through the control bottle cap.

### QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

### LIMITATIONS OF PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator. For professional use only.

### TECHNICAL ASSISTANCE

**In the USA**, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

**Outside of the USA**, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
AAC-MP 286-606	Alcohol/Ammonia, Multi-Pack Dropper Tips	6 x 3.5 mL Pkg. 100

**LOT** AA24011A  
AA24012A

 2024-01-31

**thermo**  
scientific

For Insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

# MAS<sup>®</sup> Alcohol/Ammonia

## LIQUID ASSAYED ALCOHOL/AMMONIA CONTROL

**IVD**

	<b>CON</b>						<b>SI</b>					
	AA24011A		AA24012A		<b>U</b>	AA24011A		AA24012A		<b>U</b>		
	$\bar{x}$	<R>	$\bar{x}$	<R>		$\bar{x}$	<R>	$\bar{x}$	<R>			
<b>AMMONIA</b>												
Abbott Aeroset/Architect Systems Glutamate Dehydrogenase (6K89)	69.7	55.8 - 83.7	267	213 - 320	µg/dL	41.0	32.8 - 49.1	157	125 - 188	µmol/L		
Abbott Alinity Systems Glutamate Dehydrogenase (08P22)	70.5	56.4 - 84.6	258	206 - 309	µg/dL	41.4	33.1 - 49.7	151	121 - 182	µmol/L		
Beckman Coulter AU Series Glutamate Dehydrogenase (OSR61154)	87.0	57.9 - 116	338	270 - 405	µg/dL	51.1	34.0 - 68.2	198	159 - 238	µmol/L		
Beckman Coulter AU Series Glutamate Dehydrogenase (TR60101)	**		**		µg/dL	**		**		µmol/L		
Beckman Coulter DXC/Synchron Systems Glutamate Dehydrogenase (AMM)	**		**		µg/dL	**		**		µmol/L		
Ortho-Clinical Vitros Chemistry Systems Bromophenol Blue (AMON) <sup>(a)</sup>	**		**		µg/dL	**		**		µmol/L		
Roche Cobas c/Modular Systems Glutamate Dehydrogenase (NH3L / NH3L2)	86.0	68.8 - 103	328	262 - 393	µg/dL	50.5	40.4 - 60.6	193	154 - 231	µmol/L		
Roche Hitachi/Modular D & P Systems Glutamate Dehydrogenase	**		**		µg/dL	**		**		µmol/L		
Roche Integra Systems Glutamate Dehydrogenase	86.9	69.5 - 104	325	260 - 390	µg/dL	51.0	40.8 - 61.2	191	153 - 229	µmol/L		
Siemens Advia Chemistry Systems Glutamate Dehydrogenase (AMM)	**		**		µg/dL	**		**		µmol/L		
Siemens Atellica Systems Glutamate Dehydrogenase (Amm)	93.5	60.5 - 127	366	293 - 439	µg/dL	54.9	35.5 - 74.3	215	172 - 258	µmol/L		
Siemens Dimension Systems Glutamate Dehydrogenase (AMM)	**		**		µg/dL	**		**		µmol/L		
Siemens Dimension Vista Glutamate Dehydrogenase (AMM)	96.7	77.4 - 116	349	279 - 419	µg/dL	56.8	45.4 - 68.2	205	164 - 246	µmol/L		

# MAS<sup>®</sup> Alcohol/Ammonia

## LIQUID ASSAYED ALCOHOL/AMMONIA CONTROL

**IVD**

	<b>CON</b>						<b>SI</b>					
	AA24011A		AA24012A		<b>U</b>	AA24011A		AA24012A		<b>U</b>		
	$\bar{x}$	<R>	$\bar{x}$	<R>		$\bar{x}$	<R>	$\bar{x}$	<R>			
<b>ETHANOL</b>												
Abbott Aeroset/Architect Systems												
Alcohol Dehydrogenase (3L36)	39.6	31.7 - 47.5	145	116 - 174	mg/dL	8.59	6.87 - 10.3	31.4	25.2 - 37.7	mmol/L		
Abbott Alinity Systems												
Alcohol Dehydrogenase (08P41)	40.2	32.2 - 48.3	151	121 - 181	mg/dL	8.73	6.99 - 10.5	32.7	26.2 - 39.3	mmol/L		
Beckman Coulter AU Series												
Alcohol Dehydrogenase (DiaSys)	**		**		mg/dL	**		**		mmol/L		
Beckman Coulter AU Series												
Alcohol Dehydrogenase (OSR9K229)	43.2	34.6 - 51.8	154	124 - 185	mg/dL	9.38	7.50 - 11.3	33.5	26.8 - 40.2	mmol/L		
Beckman Coulter DXC/Synchron Systems												
Alcohol Dehydrogenase (ETOH)	**		**		mg/dL	**		**		mmol/L		
Ortho-Clinical Vitros Chemistry Systems												
Alcohol Dehydrogenase (ALC)	**		**		mg/dL	**		**		mmol/L		
Roche Cobas c/Modular Systems												
Alcohol Dehydrogenase (ETOH2)	42.5	34.0 - 51.0	151	121 - 181	mg/dL	9.22	7.38 - 11.1	32.8	26.3 - 39.4	mmol/L		
Roche Hitachi/Modular D & P Systems												
Alcohol Dehydrogenase (DRI)	**		**		mg/dL	**		**		mmol/L		
Roche Hitachi/Modular D & P Systems												
Alcohol Dehydrogenase (Roche)	**		**		mg/dL	**		**		mmol/L		
Roche Integra Systems												
Alcohol Dehydrogenase	35.9	28.7 - 43.1	130	104 - 156	mg/dL	7.79	6.24 - 9.35	28.3	22.7 - 34.0	mmol/L		
Siemens Advia Chemistry Systems												
Alcohol Dehydrogenase (ETOH_2)	**		**		mg/dL	**		**		mmol/L		
Siemens Atellica Systems												
Alcohol Dehydrogenase (ETOH)	43.8	34.9 - 52.7	156	125 - 187	mg/dL	9.51	7.58 - 11.4	33.9	27.1 - 40.7	mmol/L		
Siemens Dimension Systems												
Alcohol Dehydrogenase (ETOH)	40.2	32.2 - 48.2	144	116 - 173	mg/dL	8.73	6.98 - 10.5	31.4	25.1 - 37.6	mmol/L		
Siemens Dimension Vista												
Alcohol Dehydrogenase/UV (ETOH)	42.2	33.8 - 50.6	155	124 - 186	mg/dL	9.16	7.33 - 11.0	33.7	27.0 - 40.5	mmol/L		
Thermo Fisher Scientific Indiko Systems												
Alcohol Dehydrogenase (DRI)	**		**		mg/dL	**		**		mmol/L		

**Architect®**, **Aeroset®**, Reg. TM: Abbott Laboratories, Inc., Abbott Park, IL  
**Advia® Systems**, Reg. TM: Siemens Healthcare Diagnostics, Tarrytown, NY  
**AU Series®**, **Synchron®/UniCel DxC® Series**, Reg. TM: Beckman-Coulter Instruments, Inc., Brea, CA  
**Dimension®**, **Dimension Vista®**, Reg. TM: Siemens Healthcare Diagnostics, Glasgow, DE  
**Vitros®**, Reg. TM: Ortho-Clinical Diagnostics, Rochester, NY  
**Cobas®**, **Hitachi®**, **Modular®**, **Integra®**, Reg. TM: Roche Diagnostics, Indianapolis, IN

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

**If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.**

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**USA**  
☎ 800-232-3342  
sales.diagnostics.fmt@thermofisher.com

**China**  
☎ +86 800-810-5118  
cdx.cn.info@thermofisher.com

**New Zealand**  
☎ 0800 933 966  
auinfo@thermofisher.com

**Switzerland & Austria**  
☎ +41 26 663 86 70  
cdx.ch.info@thermofisher.com

**Asia Pacific**  
☎ +61 1800 333 110  
cdd.asia.info@thermofisher.com

**France**  
☎ +33 1 40 86 65 20  
cdx.fr.info@thermofisher.com

**Nordic**  
☎ +47 2 325 0433  
info.nordic.cdd@thermofisher.com

**United Kingdom & Ireland**  
☎ +44 1442 868 940  
cdx.uk.info@thermofisher.com

**Australia**  
☎ +61 1800 333 110  
auinfo@thermofisher.com

**Germany**  
☎ +49 0800-40 40 771  
cdx.de.info@thermofisher.com

**South Africa**  
☎ +27117926790  
support-za.idd@thermofisher.com

**Canada**  
☎ 800-282-4075  
info.cddcanada@thermofisher.com

**Japan**  
☎ +81 (0)120-147-075  
JPYOK-CDD.QC@thermofisher.com

**Spain, Portugal & Italy**  
☎ +34 93589 8338  
cdx.es.info@thermofisher.com

**For countries not listed:**  
distributor.cdd@thermofisher.com

**Microgenics Corporation**  
46500 Kato Road  
Fremont, CA 94538 USA  
ISO 13485 Certified Company  
www.thermofisher.com/diagnostics

**EC REP**

**B-R-A-H-M-S GmbH**  
Neuendorfstrasse 25  
16761 Hennigsdorf, Germany  
Tel: +49 (0) 800 404 077 12  
Fax: +49 (0) 800 404 077 13

AAC-INS-VA  
Rev. 17 2020 12

**thermo**  
scientific