

# FECAL OCCULT BLOOD (FOB) QUALITY CONTROL



## NON ASSAY SPECIFIC CONTROLS INTENDED FOR PROFESSIONAL USE

Highly customised control targeting low positive levels to ensure the sensitivity of your method as well as to monitor the precision of your laboratory method and technique.

## THIRD PARTY QUALITY CONTROL

Intended to determine functionality and sensitivity of your Faecal Occult Blood test cassette.



LYOPHILIZED  
HUMAN HAEMO-  
GLOBIN



LONGER  
SHELF LIFE



EASY  
OPERATION

# PRODUCT INFORMATION

ITEM	QUANTITY	STORAGE & STABILITY
ASE FOB Control (0.5ml)	12 vials/ kit	<b>BEFORE RECONSTITUTION:</b> 2-8°C    36 months  <b>AFTER RECONSTITUTION:</b> 2-8°C    7 days

## PREPARATION GUIDELINE

1



**Preparation:**  
 Reconstitute the QC material with exactly 0.5ml deionised water.



### COMPONENTS


Lyophilized FOB Control	12 vials
Package Insert	1 pc

2



**Mix QC material and reagent:**  
 If FOB cassette diluent is 2ml  
**0.5ml QC + 2ml reagent**  
 If FOB cassette diluent is 1ml  
**0.25ml QC + 1ml reagent**

3



Drop the mixed solution onto the cassette well to view results.

## RESULT CONCENTRATION

Mixing ASE Fecal Occult Blood (FOB) Quality Control (QC) with the FOB diluent buffer in the test kit will result in a concentration of:



**15-25ng/ml**  
 (below cut off value)




**75-100ng/ml**  
 (above cut off value)



**245-265ng/ml**

\*concentration for 0.5ml FOB QC + 2ml FOB Diluent Buffer

## ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY	<b>MDA APPROVED</b>  Certificate No.: IVDB51555384919
94.36.5	ASE FOB Positive Control	1 kit	
94.16.5	ASE FOB Negative Control	1 kit	

### DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

**Address:** B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.  
**Tel:** +603-5021 4600 / 4677    **Fax:** +603-5021 4668    **E-mail:** inquiry@daiichibiotech.com.my

