

FULL BLOOD COUNT (FBC) QUALITY CONTROL

R&D SYSTEMS™
CLINICAL CONTROLS




WHOLE BLOOD QUALITY CONTROL FOR CLINICAL LABORATORY USE

A product you can rely on, R&D Systems ensure your laboratory testing routine goes smoothly. Produced over 45 different products for all major hematology instruments manufactured throughout the world, this product line also includes:

- controls
- calibrators
- linearity
- survey materials

THREE PARTS & FIVE PARTS HEMATOLOGY QUALITY CONTROL

Tri-level control for monitoring three part and five part differential instruments as well as manual methods. Each quality control is specifically designed to be compatible with every analyzer models to suit each and every needs of the laboratory use.










THREE LEVELS QUALITY CONTROL

Each kit contains 6 vials of LOW, NORMAL & HIGH level Quality Control (2 vials each).









Disks are available for uploading assay values for more efficient operation.

PRODUCT INFORMATION

ITEM	VOLUME	STORAGE	OPEN VIAL	CLOSED VIAL
5 PART QUALITY CONTROL				
 CBC X	6 x 4.5ml	2°C – 8°C	15 days	75 days
 CBC 3K	6 x 3.0ml	2°C – 8°C	8 days	75 days
 CBC 5D	6 x 4.0ml	2°C – 8°C	14 days	105 days
 CBC 5DMR	6 x 3.0ml	2°C – 8°C	14 days	75 days
3 PART QUALITY CONTROL				
 CBC ST PLUS	6 x 2.5ml	2°C – 8°C	14 days	105 days
 CBC TECH	6 x 3.5ml	2°C – 8°C	14 days	75 days
 CBC 3D	6 x 3.0ml	2°C – 8°C	14 days	105 days

COMPATIBLE ANALYZERS

ABBOT CELL-DYN

 RUBY  EMERALD  SAPPHIRE  1700  1800  3200  3500  3700



BECKMAN COULTER

 DXH 800  HMX  GEN-S  LH 500  LH 700

HORIBA

 ABX Micros ES60





MINDRAY

 BC 6600  BC 6800  BC 5000 series  BC 3000 series  BC 30s





NIHON KOHDEN

 CELLTAC 6000  CELLTAC 7000  CELLTAC 8000

SIEMENS

 ADVIA 120  ADVIA 2120  ADVIA 2120i  ADVIA 60

SYSMEX

 XN series  XT series  XE series  XS-800  KX-21  XP-100

ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY
X003X	CBC X	3 Levels / kit
3K303X	CBC 3K	3 Levels / kit
5D004	CBC 5D	3 Levels / kit
5DMR02	CBC 5DMR	3 Levels / kit
ST208	CBC ST Plus	3 Levels / kit
CT003X	CBC Tech	3 Levels / kit
3D506	CBC 3D	3 Levels / kit



DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

Address: B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.
Tel: +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my

