

# ERYTHROCYTE SEDIMENTATION RATE (ESR) QUALITY CONTROL

**R&D SYSTEMS™**  
CLINICAL CONTROLS



## INTENDED FOR MANUAL AND AUTOMATED ESR METHOD

SEDRite Plus is a bi-level control formulated to provide values in the clinically normal and elevated ranges. The control is designed to monitor Erythrocyte Sedimentation Rate (ESR) values obtained from manual and automated ESR methods. It is also used in laboratories to monitor the factors that cause variability in ESR results, such as:

- technique
- time
- temperature
- tube position

## CLOSE TO HUMAN SAMPLE

SEDRite Plus is an excellent control for the ESR tests because the control cells rouleaux in the same manner as fresh whole blood cells.

## INDEPENDENT ASSESSMENT

Third party whole blood count differential control to prevent result bias.



## 2 LEVELS QUALITY CONTROL

Formulated to provide values in the clinically normal and elevated ranges.

Disks are available for uploading assay values for more efficient operation.

## PRODUCT INFORMATION

ITEM	VOLUME	STORAGE	OPEN VIAL	CLOSED VIAL
SEDRITE PLUS	4 x 9.0ml	2°C – 8°C	30 days	195 days

## COMPATIBLE ANALYZERS

### BECTON DICKINSON

- SED I System Geo II V.3xx
- SED I 15

### BIOCADE

- Bioca VS12/40

### BIOCODE HYCEL

- Lena
- Eriline
- Therna NE
- Biostatic

### CLINICAL DATA

- Excyte 10/M
- Excyte 40

### COFIM

- SR50/60

### DIESSE

- Mini-Ves
- Vesmatic 10/20/30
- Vesmatic Cube 30/80/200
- Cube 30 Touch

### ELECTA LAB

- MicroSed
- Monitor 20/20E/100

### ELITECH

- Mixrate 10/20/100

### GREINER/SAYAG

- SRS-20/100
- SRT/SRTX

### HORIBA

- Auto Compact Starrsed
- Starrsed Inversa 24M

### KIMASED

- Auto 16/20/60

### MANUAL METHOD

- WESTERGREN  
saline diluted  
sodium citrate diluted  
undiluted
- WINTROBE

### OTHERS

- LS DIAG SRA 48LS
- ALL DIAG VES STATIC
- BERKHUN SDM-60
- Other instruments

### ROCHE

- InteRRliner
- Starliner

### STARSED I

- SediPlus S100/200/2000

### TERUMO

- Monitor V10/20/100
- Monitor J

### VGA/EUROTECH

- VT-25/50/72/100

### XUNDA

- ESR30

## ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY
SR002X	Sedrite Plus	2 Levels / kit

**MDA APPROVED**



### DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

**Address:** B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.  
**Tel:** +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my

