

CBC MONITOR ONLINE PROGRAM

INTERNAL QUALITY CONTROL (IQC) PROGRAM

MONITOR YOUR OWN PEER GROUP PERFORMANCE

The CBC Monitor Online Program is available free of charge to all users. Simply upload or update the data in the program and peer group report can be generated. Disks are available for uploading assay values to analyzers.

CBC-Monitor2 includes the features of the CBC-Monitor program plus the ability to track comparison history for a full year instead of 6 months. The program is available in English and French. Customers may enroll in the program anytime, without making any phone calls.

**SIGN UP
FREE!**

✓ PEER GROUP REPORT


Get your peer report from a group of similar instrument types, using the same control product of the same lot and level

✓ REMOTE MONITORING

View your data reports by monthly, yearly and among peers- anytime, anywhere

✓ VIEW RESULTS YOUR OWN WAY

Choose statistical results to view Peer Comparison, Quality Index Graphs, Levey-Jennings Charts, Peer History, and compare with other instrument



1 Upload/fill in data values

CBC-Monitor
TerumoV10, V20, V100, Tec, Hycl
SR075

Peer History


Year	2015		2015		2015		2015		2015		2015		2014	
	Aug	Jul	Jun	May	Apr	Mar	Feb	Jan	Dec	Nov	Oct	SR075	SR074	
ESR	24	31	2	2	2	2	29	29						
Level 1														
Number of Labs	3.8	3.0	3.7	3.0	3.5	4.3	4.0	4.0						
Your Lab Mean	4.3	4.3	3.7	3.0	3.9	4.3								
Group Mean	-0.5	-1.3	0.0	0.0	-0.4									
Mean Difference	1.2	1.2	1.5	1.4	1.4									
Your Lab SD	1.5	1.6	1.6	1.6	1.6	1.4								
Group SD	-0.34	-0.82	0.02	0.00	-0.27									
SDI	31.4	40.0	40.3	46.7	39.6									
Your Lab %CV	35.0	37.4	42.5	51.7	41.5	33.3								
Group %CV	0.90	1.07	0.95	0.90	0.95									
PI														

CBC-Monitor
SR075

Other Instruments

Peer Group	Level 1			Level 2		
	Number of Labs	Group Mean	Group SD	Number of Labs	Group Mean	Group SD
Biocad, Electalab, TerumoV10, V20, V100, Tec, Hycl, ESR	24	4.3	1.5	35.0	24	62.8
					9.0	14.3

Peer Group	Level 1			Level 2		
	Number of Labs	Group Mean	Group SD	Number of Labs	Group Mean	Group SD
EUROTEC Med Syst VT-72 - UNDILUTED - World	11	5.9	1.5	24.6	11	26.2
					9.3	35.4



2 Generate reports for monitoring