

BIOCHEMISTRY QUALITY CONTROL

ThermoFisher
SCIENTIFIC

THERMO SCIENTIFIC MAS



Thermo Scientific MAS Quality Control – your ultimate solution for precision and accuracy in laboratory testing. Designed with cutting-edge technology and trusted by professionals worldwide, this quality control system ensures consistent and reliable results every time. With its user-friendly interface and comprehensive data analysis capabilities, you can streamline your workflow and make informed decisions with ease. Elevate your laboratory standards and unleash the true potential of your research.



SERUM CHEMISTRY CONTROLS

- MAS ChemTRAK•H
- MAS Alcohol/Ammonia
- MAS Bilirubin

URINE CHEMISTRY CONTROLS

- MAS UA Control
- MAS UrichemTRAK

TOXICOLOGY CHEMISTRY CONTROLS

- MAS DOA TOTAL
- MAS TOX Control

IMMUNOASSAY CONTROLS

- MAS Liquimmune
- MAS PAR TDM
- MAS T-Marker

IMMUNOLOGY & PROTEIN CONTROLS

- MAS Immunology
- MAS CSF

CARDIAC & DIABETES CONTROLS

- MAS CardiImmune-XL
- MAS Diabetes

	TYPE OF TESTS	QUALITY CONTROL	CAT. NO.	LEVEL	SIZE	STABILITY
O R D E R I N G I N F O R M A T I O N S	SERUM CHEMISTRY QUALITY CONTROL					
	General Chemistry /TDM <i>human serum</i>	MAS chemTRAK -H <i>(Assayed)</i>	CHA-1 CHA-2 CHA-3	Level 1 Level 2 Level 3	6 x 5ml 6 x 5ml 6 x 5ml	Closed vial 30 months -25°C to -15°C
		MAS chemTRAK -H <i>(Unassayed)</i>	CHU-1 CHU-2 CHU-3	Level 1 Level 2 Level 3	10 x 15ml 10 x 15ml 10 x 15ml	Open vial 7 days 14 days 2-8°C (<i>Assayed</i>) 2-8°C (<i>Unassayed</i>)
Ethanol/Ammonia <i>human serum</i>	MAS Alcohol/Ammonia	AAC-MP	Bi level (<i>Multipack</i>)	6 x 3.5ml	Closed vial 24 months 2-8°C Open vial 30 days 2-8°C	
Bilirubin Test <i>bovine serum</i>	MAS Bilirubin	BC101 BC102 BC103	Level 1 Level 2 Level 3	6 x 5ml 6 x 5ml 6 x 5ml	Closed vial 24 months 2-8°C Open vial 14 days 2-8°C	
URINE CHEMISTRY QUALITY CONTROL						
Urine FEME <i>human urine</i>	MAS UA	MAUAB-115 MAUAB-215 MAUAB-160 MAUAB-260 MAUAB-MP	Level 1 (<i>Normal</i>) Level 2 (<i>Abnormal</i>) Level 1 (<i>Normal</i>) Level 2 (<i>Abnormal</i>) Bi level (<i>Multipack</i>)	6 x 15ml 6 x 15ml 4 x 60ml 4 x 60ml 6 x 15ml	Closed vial 24 months 2-8°C Open vial 3 months 6 weeks 2-8°C 18-25°C	
Urine Chemistry <i>human urine</i>	MAS UrichemTRAK	UR11001 UR22002 UR-MP	Level 1 Level 2 Bi level (<i>Multipack</i>)	6 x 15ml 6 x 15ml 6 x 15ml	Closed vial 24 months 2-8°C Open vial 24 months 2-8°C	
IMMUNOASSAY QUALITY CONTROL						
Immunoassay Test <i>human serum</i>	MAS Liquimmune	LIG-101 LIG-102 LIG-103	Level 1 Level 2 Level 3	6 x 5ml 6 x 5ml 6 x 5ml	Closed vial 60 months 90 days -25°C to -15°C 2-8°C Open vial 30 days 2-8°C	
Therapeutic Drugs /Immunoassay <i>bovine serum</i>	MAS PAR TDM	PTD1-1001 PTD1-2002 PTD1-3003	Level 1 Level 2 Level 3	6 x 5ml 6 x 5ml 6 x 5ml	Closed vial 36 months 2-8°C Open vial 30 days 2-8°C	
Cancer Antigen Tests <i>human serum</i>	MAS T marker	TUM-101 TUM-202 TUM-303	Level 1 Level 2 Level 3	6 x 3ml 6 x 3ml 6 x 3ml	Closed vial 36 months 90 days -25°C to -15°C 2-8°C Open vial 30 days 2-8°C	
CARDIAC & DIABETES QUALITY CONTROL						
Serum Cardiac Marker <i>human serum</i>	MAS Cardioimmune XL	CAI-XLL CAI-XL1 CAI-XL2 CAI-XL3 CAI-XL4	Level L Level 1 Level 2 Level 3 Tri level (<i>Multipack</i>)	6 x 3ml 6 x 3ml 6 x 3ml 6 x 3ml 6 x 3ml	Closed vial 36 months -25°C to -15°C Open vial 30 days 15 days 2-8°C 2-8°C	
Hemoglobin A1c <i>human whole blood</i>	MAS Diabetes	DBCL-MP	Bi level (<i>Multipack</i>)	6 x 1ml	Closed vial 24 months 60 days -25°C to -15°C 2-8°C Open vial 30 days 2-8°C	
IMMUNOLOGY & PROTEIN QUALITY CONTROL						
Serum Protein /Serology <i>human serum</i>	MAS Immunology	IC-101 IC-202 IC-303	Level 1 Level 2 Level 3	6 x 3ml 6 x 3ml 6 x 3ml	Closed vial 24 months 2-8°C Open vial 30 days 2-8°C	
Cerebrospinal Fluid <i>human serum</i>	MAS CSF	CSF-MP	Bi level (<i>Multipack</i>)	6 x 3.5ml	Closed vial 24 months 2-8°C Open vial 30 days 2-8°C	
TOXICOLOGY QUALITY CONTROL						
Urine Toxicology <i>human urine</i>	MAS DOA TOTAL	DOAT-1 DOAT-2 DOAT-3 DOAT-4 DOAT-5 DOAT-6 DOAT-MP	Level 1 Level 2 Level 3 Level 4 Level 5 Level 6 6 Levels (<i>Multipack</i>)	6 x 18ml 6 x 18ml 6 x 18ml 6 x 18ml 6 x 18ml 6 x 18ml 6 x 18ml	Closed vial 24 months 2-8°C Open vial 30 days 2-8°C	
Serum/Plasma Toxicology <i>bovine plasma</i>	MAS TOX	10011608	Tri level (<i>Multipack</i>)	6 x 5ml	Closed vial 24 months 2-8°C Open vial 30 days 2-8°C	