

HIV 1.2.O RAPID TEST CASSETTE

Beright[™]



FAST RESULT



EASY
TO USE



10
MINUTES
RESULT



YOUR RELIABLE COMPANION IN THE FIGHT AGAINST HIV

This innovative chromatographic immunoassay is specifically designed to detect antibodies of Human Immunodeficiency Virus (HIV) type 1, type 2, and subtype O in whole blood, serum or plasma. Its advanced technology ensures high sensitivity and specificity, minimizing the chances of false positives or negatives.


By providing early detection capabilities, this test allows for timely intervention and support. It plays a crucial role in preventing the spread of HIV by enabling prompt treatment and counseling.

SENSITIVITY
99.9 %

SPECIFICITY
99.9 %

PRODUCT INFORMATION

ITEM	QUANTITY	STORAGE & STABILITY
Beright HIV 1.2.O Rapid Test Cassette	40 Tests/ kit	2°C – 30°C shelf life: 23 months




**WHOLE BLOOD
/ SERUM
/ PLASMA**

COMPONENTS

Test Cassettes	40pcs
Droppers	40pcs
Buffer	2pcs
Package insert	1pc

RESULT INTERPRETATION

POSITIVE




Two or three distinct colored lines appear. One line should always appear in the control line region (C).

NEGATIVE



No colored line appears in the test line regions (T1 and T2).

INVALID



Control line fails to appear.

ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY
IHI-T402	Beright HIV 1.2.O Rapid Test Cassette	40 Tests / kit



DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

Address: B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.
Tel: +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my

