

HBsAg RAPID TEST CASSETTE

Beright[™]



FAST RESULT



EASY
TO USE



15
MINUTES
RESULT



SAFEGUARD YOUR PATIENTS' HEALTH AGAINST HEPATITIS B

Designed for professional in vitro diagnostic use, this rapid test provides a quick and accurate qualitative detection of Hepatitis B Surface Antigen (HBsAg) in whole blood, serum, or plasma. With its advanced chromatographic immunoassay technology, healthcare professionals can confidently diagnose and monitor patients for this viral infection.

This test cassette offers convenience and efficiency, making it an indispensable tool in healthcare settings. Its user-friendly design allows for easy handling and interpretation of results, saving precious time for medical personnel. The clear and distinct test line and control line ensure accurate readings, leaving no room for interpretation errors.

SENSITIVITY
99.3 %

SPECIFICITY
99.8 %

PRODUCT INFORMATION

ITEM	QUANTITY	STORAGE & STABILITY
Beright HBsAg Rapid Test Cassette	40 Tests/ kit	2°C – 30°C shelf life: 23 months



**WHOLE BLOOD
/ SERUM
/ PLASMA**

COMPONENTS

Test Cassettes	40pcs
Droppers	40pcs
Buffer	2pcs
Package insert	1pc

RESULT INTERPRETATION

POSITIVE



Two distinct colored lines appear. One colored line should be in the control region (C) and another colored line should be in the test region (T).

NEGATIVE



No colored line appears in the test line regions (T).

INVALID



Control line fails to appear.

ORDERING INFORMATION

CAT. NO	ITEM	QUANTITY
IHBSG-402	Beright HBsAg Rapid Test Cassette	40 Tests / kit

MDA APPROVED



Certificate No.: IVDD1493421-53277

DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

Address: B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.
Tel: +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my

