

# LC-560 & LC-570 PNEUMATIC SHOCKWAVE

## Lifotronic

LIGHT  
**900G**  
LESS  
FATIGUE

POWERFUL  
**22 HZ**  
/5 BARS

EFFICIENT  
**3.0**  
MIN/  
SESSION

**12**  
MONTHS  
WARRANTY



### NON-INVASIVE

Extracorporeal Shock Wave Therapy (ESWT) is a non-invasive therapy used to treat pain and promote healing by delivering high-energy sound waves to injured or damaged tissue during short treatment sessions.

Through the use of targeted ESWT, inflamed tissue that causes pain will begin to regenerate because of increased blood circulation created by the treatments. Many patients see an improvement in their pain levels after as few as one or two sessions of ESWT.

### FEATURES

- Real-time pressure indicating**  
Parameters can adjust during therapy
- Shock-absorbing applicator**  
Avoid injury to operators
- Spinning projectile**  
Better collimation and higher durability

# ACCESSORIES

## SHOCKWAVE TRANSMITTER



Standard Transmitter



Tendon Transmitter



Finger Transmitter



Point Transmitter



Deep Transmitter



Trigger Transmitter



Larger Transmitter

## MESSAGE HEAD



Massage gun

Spherical head

Flat type head

U shaped head

Cylindrical head

ESWT is beneficial for patients because there is no down time after the procedure, so time missing for work, school or sports practice is not an issue for patients.

## PORTABLE ACCESSORIES



Tablet

Trolley

# VERSATILE FOR MULTIPLE DISCIPLINES

## CHRONIC TENDINOPATHIES


- ☑ Calcifying tendinopathy of the shoulder
- ☑ Lateral epicondylopathy of the elbow (tennis elbow)
- ☑ Greater trochanter pain syndrome (GTPS)
- ☑ Patellar tendinopathy
- ☑ Achilles tendinopathy
- ☑ Plantar fasciitis, with or without heel spur

## BONE PATHOLOGIES

- ☑ Delayed bone healing
- ☑ Stress fracture
- ☑ Bone Non-Union (pseudarthrosis)
- ☑ Avascular bone necrosis without articular derangement
- ☑ Osteochondritis Dissecans (OCD) without articular derangement

## SKIN PATHOLOGIES

- ☑ Non-circumferential burn wounds
- ☑ Skin ulcers
- ☑ Delayed or non-healing wounds



## HAND-PIECE OPERATION

- ☑ Uninterrupted pressure adjustment avoiding unwanted turns and moves
- ☑ Maximize patient comfort
- ☑ Optimized workflow
- ☑ Maximum operator comfort



## SPECIFICATION

<b>Model</b>	LC-560: 1 handpiece, LC-570 : 2 handpieces
<b>Power</b>	Pneumatic
<b>Pressure range</b>	1.0-5.0 bar
<b>Frequency</b>	1-22 Hz
<b>Shock count range</b>	100-9900
<b>Interface</b>	LCD screen and mechanical keys
<b>Input power</b>	240 VA
<b>Controller Dimensions</b>	15cm(H) x 43cm(W) x 32cm(D)
<b>Controller Weight</b>	13.3 kg
<b>Power Supply</b>	100-240V AC, 50/60 Hz
<b>Optional accessories</b>	Tablet, massage gun, trolley

MDA APPROVED



Certificate No.: GB9425023-135687

### DAIICHI BIOTECH SERVICES SDN BHD (925272-K)

**Address:** B-SG-20, Sunway Geo Avenue, Jalan Lagoon Selatan, Bandar Sunway, 47500 Subang Jaya, Selangor.

**Tel:** +603-5021 4600 / 4677 **Fax:** +603-5021 4668 **E-mail:** inquiry@daiichibiotech.com.my



daiichibiotech.com.my